

# Western Instruments

Established 1965

26509 Township Road 543  
Sturgeon County, Alberta  
T8T 1M1

Phone: 780-459-6720  
Fax: 780-459-7837  
E-mail: [info@westerninstruments.com](mailto:info@westerninstruments.com)

Web: [www.westerninstruments.com](http://www.westerninstruments.com)

## Ultrasonic Couplant

Western Instruments Ultrasonic Couplant that has been developed for our Northern Climate and local process industries. The Couplant has been descriptively named; *-40° Ultrasonic Couplant*.

### *-40° Ultrasonic Couplant*

This Couplant is ideal for use with Portable UT Equipment year round. Obviously this material is used extensively in our Cold winter months. However, it maintains its low viscosity above +40°C / 105°F as well. Other features of our *-40° Ultrasonic Couplant* are the addition of a UV tracer, to ensure its full removal can be inspected. *-40° Ultrasonic Couplant* is removed with a dampened cloth or industrial wipe. For flaw detection, *-40° Ultrasonic Couplant* can be diluted with water to raise its viscosity.

*-40° Ultrasonic Couplant* is offered in the following packaging;

- 100ml / 4oz. Squeeze Bottle.
- 500ml / 16oz. Squeeze Bottle.
- 3.8L / 1 Gallon Kit, Complete with pump, and 2 Squeeze Bottles.
- 19L / 5 Gallon Pail.
- 205L / 45 Gallon Drum.



Ultrasonic Couplant 4oz (100ml)